KYERON

Doc: IFU02-ELIZ v4.1 (2025-06-19)

KYERON

ELIZ[®]

Synthetic Bone Substitute, Granules

Synthetic Bone Substitute, Block

EN	INSTRUCTIONS FOR USE
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
CS	NÁVOD K POUŽITÍ
DA	BRUGERVEJLEDNING
DE	GEBRAUCHSANWEISUNG
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ES	INSTRUCCIONES DE USO
ET	KASUTUSJUHEND
FI	KÄYTTÖOHJEET
FR	MODE D'EMPLOI
HR	UPUTE ZA UPORABU
HU	HASZNÁLATI UTASÍTÁS
IS	NOTKUNARLEIÐBEININGAR
IT	ISTRUZIONI PER L'USO
LT	NAUDOJIMO INSTRUKCIJA
LV	LIETOŠANAS INSTRUKCIJAS
NL	GEbruIKSAANWIJZING
NO	BRUKSANVISNING
PL	INSTRUKCJE STOSOWANIA
PT	INSTRUÇÕES DE UTILIZAÇÃO
RO	INSTRUCȚIUNI DE UTILIZARE
SK	NÁVOD NA POUŽITIE
SL	NAVODILA ZA UPORABO
SV	BRUKSANVISNING

Ceramisis Limited,
914 Herries Road, Sheffield, S6 1QW, UK.
e-mail: info@ceramisis.com

Medical Device Management Ltd
Block B, The Crescent Building, Northwood,
Santry, Dublin, D09 C6X8, Ireland

EN

INSTRUCTIONS FOR USE ELIZ®

Synthetic bone substitute, granules
Synthetic bone substitute, block

Material & Performance Characteristics

ELIZ® is a resorbable porous material in block and granular form that is composed of nominally 60% (w/w) hydroxyapatite and 40% (w/w) β -tricalcium phosphate. The device will perform as an osteoconductive scaffold for at least 12 months, following which the residual material will continue to slowly resorb mainly within a period of 24 months from implantation. Residual fragments may remain that will continue to be slowly resorbed over time as part of the patient's natural bone remodelling.

Intended Purpose

ELIZ® is a synthetic bone substitute intended to fill traumatically or surgically created bone defects that are not intrinsic to the stability of the bony structure.

Indications-for-use

ELIZ® is indicated for treating defects of the long bones and extremities. ELIZ® granules may also be used for interbody fusion of the spine.

Clinical Benefit

ELIZ® supports the ingrowth and fusion of adjacent viable bone when placed within an osseous defect and provides a synthetic alternative or supplement to human or animal cadaver bone that is free from the risk of disease transmission and can in many cases spare the patient the trauma of autograft harvesting.

Contraindications

ELIZ® is contraindicated where the device is intended to provide structural support in the skeletal system and **MUST NOT BE USED TO GAIN SCREW FIXATION**. ELIZ® is also contraindicated in patients with:

- Acute or chronic infections, particularly at the implantation site
- Severe vascular or neurological disease
- Uncontrolled diabetes
- Severe degenerative disease
- Collagen disease
- Hypercalcaemia, abnormal calcium metabolism
- Inflammatory bone disease
- Malignant tumours
- Severely impaired renal function
- Defects of open epiphyseal plates

Warnings

The device is provided **STERILE** by gamma irradiation, and within sealed double packaging to allow aseptic presentation. **DO NOT USE IF OPENED, PUNCTURED OR PRODUCT APPEARS DAMAGED**. Read the expiration date before use and **DO NOT USE BEYOND THE EXPIRATION DATE**.

ELIZ® is for **SINGLE USE ONLY**. Do not attempt to re-sterilise or re-use due to the risk of cross infection. Do not cut or reshape the implant whilst implanted in the surgical site to prevent particulate debris being introduced into the surgical site. Blocks & granules must be secured to prevent any potential migration of the implant and should only be used in procedures where the implant can be adequately contained within the defect. ELIZ® is opaque to x-rays and may hide areas under or above the implant on a radiograph. The device is MRI safe.

Precautions

It is important to ensure that the area around the implantation site be secured mechanically with rigid fixation to provide structural support and maintain the implant in a static, load free environment. In order to facilitate the formation of new bone, ELIZ® should only be implanted in direct contact with well vascularised viable bony tissue.

The interactions with other substances are unknown but no specific examples of negative interactions have been identified.

Intended User & Intended Patient

ELIZ® is for professional use only, and is only intended for use by surgeons familiar with and skilled in the techniques of bone repair and replacement. ELIZ® is intended for use in any adult patient requiring bone graft as indicated.

Possible Complications & Adverse Reactions

Possible post-operative complications are those that may occur with any surgery, including delayed wound healing, nerve pain, neurological deficit, infection, haematoma, and pain at the implant site. A successful result may not be achieved in every case. A secondary operation to remove or replace an implant may be necessary due to surgical error, specific medical conditions, non-union, delayed or incomplete union, fracture of the newly formed bone, or device displacement fracture or crushing with or without generation of particulate debris. Allergic reaction to the product may occur in very rare circumstances.

Application

Blocks may be gently shaped with a scalpel or standard surgical bone tools to adequately fit the defect and after shaping must be rinsed in sterile water or saline to remove particulate debris before attempting to implant. Prior to implantation the implant may be soaked with blood and/or combined with autologous bone as per normal surgical practice. The device should be gently and carefully placed into the defect which should be completely filled. Secure the surgical site after implanting to prevent motion and implant displacement. If the material does not fill the site correctly, remove the implant and repeat the process with a fresh amount of ELIZ®

Storage

Store ELIZ® between 5°C and 30°C. Direct contact with sunlight or heating systems should be avoided.

Shelf Life and Disposal

The expiration date is printed on the labelling. Do not use ELIZ® after the expiration date. Residual material should be disposed as standard clinical waste.

Implant Card, Traceability & Event Reporting

The device is identifiable by its lot number which is pre-printed on the labelling and implant card. The implant card enclosed with this device should be supplied to the patient after implantation. Before providing it to the patient, the healthcare institution must complete the implant card with the following information:

 ? Name of the patient  Date of implantation

 Healthcare institution (name and address)

The summary of safety & clinical performance may be obtained via the Eudamed database <https://ec.europa.eu/tools/eudamed> where it is linked to the Basic UDI-DI: 506017162HATCPZ9, or upon request to the manufacturer.

EU: Serious incidents should be reported to the manufacturer and to the competent authority of the member state of the user and/or patient's establishment.

Other Territories: Serious incidents should be reported to the manufacturer and to the relevant national authority

Note: ELIZ® is for professional use only. Responsibility for proper selection of patients, for adequate training, for experience in the choice of ELIZ®, all aspects of the surgery, and for the choice of post-operative procedures rests entirely with the surgeon.

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