ELIZ[®]

4911025 / 4911050 / 4911100 / 4911150 / 4911200 / 4911250 / 4911300 /
4912410 / 4912412 / 4912413 / 4912414 / 4912420 / 4912430 / 4912440 / 4912450 / 4912460 / 4912470 / 4912480 / 4912490 / 4912492

Ceramisys Limited,
914 Herries Road, Sheffield S6 1QW, UK.

www.ceramisys.com/patientinfo e-mail: info@ceramisys.com

EN

IT

BG

LT

CS

LV

DA

NL

DE

NO

EL

PL

ES

PT

ET

RO

FI

RU

FR

SK

HR

SL

HU

SV

IS

EN

The product is a bone graft substitute. This means that it fills gaps in bones and supports them to heal by acting as a scaffold. The product can be used in any adult patient. The product is for use in bones of the limbs, hands, and feet. It is also for use in spine surgery (fusion of the spine).

The product contains hydroxyapatite (60%), and beta tricalcium phosphate (40%).

The product will work for at least 12 months allowing time for the bone to heal. The product will be slowly absorbed by the body, mainly within 2 years. Residual fragments may remain that will be absorbed over time. It does not need any maintenance.

Complications are those that can happen with any surgical procedure. These include delayed wound healing, nerve pain, nerve deficit, infection, hematoma and pain at the implant site. A good result may not happen in every case. Another operation may be needed. This can be caused by undiagnosed medical conditions or surgical error. It might also be needed where your bone does not heal over time. Movement or physical damage to the product after surgery might also need a second operation. Allergic reaction to the product may occur very rarely. There are no known problems of the product with other substances.

There is no special product care that you need to observe. Always follow the advice of your doctor.

The product is safe for use with magnetic resonance imaging (MRI).

Caution: The product partly blocks x-rays and can mask the area on x-ray images.

Consult your doctor if you think you have any problems after your surgery. Any serious incident that happens with the product should be reported. It should be reported to the manufacturer. It should also be reported to the health authority in your country of residence:

E.U.: to the Competent Authority of your country

Australia: to the Therapeutic Goods Administration (www.tga.gov.au)

Other Territories: to the national authority of your country

The symbols on your implant card and in this leaflet have these meanings:

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